

- Enables a fuller family life without extensive reconstruction & adaptations
- Ensures optimum social inclusion at home, school, holiday and out and about
- Safe & proven for carer & passenger alike
- Helps to overcome respite issues
- Assessment - *Free*
- Demo - *Free*
- Training - *Free*

Please call MHS: 01553 811977



Manual Handling Solutions

58, Paige Close, The Meadows, Watlington
 King's Lynn, Norfolk, PE33 0TQ
 Tel: +44 (0)1553 811977 Fax: +44 (0)1553 811004
 Mobile: +44 (0)7768 506297
 E-mail: sales@manualhandlingsolutions.co.uk
 www.mobility-mhs.co.uk
 www.manualhandlingsolutions.com



The stairclimber people

Paediatric Stairway Solutions

Maximise mobility for the whole family



www.manualhandlingsolutions.com

Full, quick release interface system compatible with most recline and tilt in space systems.

Full interface with most systems to promote full and active participation with family & friends at home, school, work or out & about by enabling easy switching between wheelchair, buggy and stairclimber.



Technical data	c-max U1	s-max
Max. lifting capacity	140kg / 160kg	120kg / 160kg
Batteries	2 x 12V 5Ah	2 x 12V 3.3Ah
DC-motor	24V	24V
Climbing speed	8-23 steps/minute (continually adjustable)	8-23 steps/minute (continually adjustable)
Capacity with one battery charge	300-500 steps (depending on the person's weight)	up to 300 steps (depending on the person's weight)
Steps height	maximum 21cm	maximum 22cm
Weight	31.7kg 17.2kg - heaviest part disassembled	22kg 14.5kg (without batteries)
Height	1090mm	1200mm
Width	485mm (total) 440mm (without arm rest)	520mm (with handle) 330mm (climbing unit)
Depth	915mm (total) 750mm (folded in foot rest)	330mm
Options available	12V/24V charging cable for the car, hip belt, seat belt system	

The AAT s-max and AAT c-max Paediatric Versions

- More mobility and quality of life. These push wheel chairs with stair climbing function allow you to overcome virtually all obstacles at home and on the road.
- The care attendant needs little physical power as these devices can be operated effortlessly and safely. People with reduced mobility are transported safely and comfortably up and downstairs.
- On stairs as well as on level ground they prove their versatility. Removable arm rests make easy shifting from one chair to another possible. Due to their compact dimensions and foldable foot rests, they are easy to manoeuvre even on very narrow stair cases or winding stairs.
- All kinds of floor coverings such as wood, stone, metal, carpet or PVC flooring can be driven on without problems.
- The device is easily and swiftly dismantled and transported in virtually any car - ready for use at your destination.
- Ideal in respite circumstances.



Bespoke units can be designed to suit customers own requirements